

Technical data

ProSim SPOT Light

Featuring an exclusive ergonomic design, the Pro-Sim SPOT Light is the first comprehensive SpO₂ functional tester to come in a handheld easy-to-use device. SPOT Light is lightweight and flexible, with three custom presets specially-designed to make it the fastest and easiest-to-use device on the market today for pulse oximeter functional testing. A helpful LCD display and three simple buttons make it effortless to rapidly change parameters and view each signal output sent to the pulse oximeter at a glance. An interchangeable, long-life battery ensures uninterrupted all-day operation without the need to connect to a power supply.

SPOT Light sets up in seconds to send SpO₂ saturation, heart rate, perfusion, transmission, artifact noise, and eight different manufacturers custom R-curves to a pulse oximeter or patient monitor. It is designed to take the abuse of everyday use, including the stress of transportation, moving from room to room and device to device, and the occasional fall from the workbench.

An increasing number of medical devices now have pulse oximetry built in. SPOT Light can be paired with other Fluke Biomedical test tools to make modern device testing task fast and easy. Need to check a defibrillator with pulse oximetry? SPOT Light with the Impulse 7000 Defibrillator Analyzer allow you to test any external defibrillator on the market today. Do you have a fleet of new patient monitors coming that require incoming inspection? Using SPOT Light with ProSim 4 Vital Signs Simulator enables you to easily accomplish this task in hours instead of days.



Key features

- Small, portable, and lightweight, weighing in at 8.8 oz
- Large, easy to read LCD display
- Rechargeable battery lasts 10 hours minimum
- Signal strength indicator
- SpO₂ Saturation 70% to 100%
- Heart rate: 30 BPM to 240 BPM
- Perfusion: 0.2%, 2% and 10%
- Transmission: dark/thick, normal and light/thin
- Artifacts: respiration and ambient light
- Three editable test sequences
- R-curves: 8 including Masimo, Nellcor and Nonin
- Expandable pairing capabilities

Specifications

| Physical specifications | |
|------------------------------|---------------------------------------------------------|
| Display | 2¼" x 1½" LCD display |
| Size (WxDxH) | 12.2 cm x 9.7 cm x 4.8 cm (4.8 in x 3.8 in x 1.9 in) |
| Weight | 8.8 oz. |
| Electrical specifications | |
| Battery | |
| Type | Rechargeable Lithium-Ion |
| Charge time | Approximately two hours |
| Operating duration | Approximately ten hours minimum |
| AC adapter | |
| Input voltage | 100 V ac to 240 V ac |
| Input frequency | 50/60 Hz |
| Input current | 0.5 A (rms) |
| Output voltage | 6 V dc |
| Environmental specifications | |
| Operating temperature | 0° C to 35° C (32° F to 95° F) |
| Storage temperature | -35° C to 50° C (-31° F to 122° F) |
| Operating humidity | 20% to 80% RH (non-condensing) |
| Measurements | |
| SpO ₂ saturation | 70%, 75%, 80%, 85%, 90%, 95%, 97%, 98%, 99%, and 100% |
| Heart Rate | 30, 40, 60, 65, 80, 100, 120, 140, 150, 180 and 240 BPM |
| Perfusion rate | 0.2%, 2% and 10% |
| Transmission level (finger) | Dark/thick, normal and light/thin |
| Artifacts | Respiration or ambient light (50 Hz or 60 Hz) |



Ordering information

Models/descriptions

| | |
|---------|-----------------------------------------------------------|
| 4111101 | ProSim SPOT Light SPO2 Functional Tester, USA |
| 4142735 | ProSim SPOT Light SPO2 Functional Tester, Schuko |
| 4142747 | ProSim SPOT Light SPO2 Functional Tester, UK |
| 4142758 | ProSim SPOT Light SPO2 Functional Tester, Japan |
| 4142764 | ProSim SPOT Light SPO2 Functional Tester, Australia/China |
| 4142814 | ProSim SPOT Light SPO2 Functional Tester, Brazil |

Standard accessories

| | |
|---------|-----------------------------------------------------------------------------------------|
| 4151274 | ProSim SPOT Light Users Manual |
| 4026823 | ProSim SPOT Light Battery Pack |
| 4026773 | ProSim SPOT Light Power Supply Line Cord ProSim SPOT Light Line Cord (country-specific) |
| 4180099 | CCPS4 SPOT Light, PS4 SPOT Light Carry Case |

Optional accessories

| | |
|---------|--------------------------------|
| 4034393 | USB Cable |
| 4026823 | ProSim SPOT Light Battery Pack |



About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance. Highly credentialed and equipped with a NVLAP Lab Code 200566-0 accredited laboratory, Fluke Biomedical also offers the best in quality and customer service for all your equipment calibration needs.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical regulatory commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 and ISO 13485 medical device certified and our products are:

- CE Certified, where required
- NIST Traceable and Calibrated
- UL, CSA, ETL Certified, where required
- NRC Compliant, where required

Fluke Biomedical

We empower our everyday heroes to focus only on protecting lives.

Fluke Biomedical
6920 Seaway Boulevard
Everett, WA 98203

For more information, contact us at:
(800) 850-4608 or Fax (440) 349-2307
sales@flukebiomedical.com
flukebiomedical.com

©2023 Fluke Biomedical.
Specifications subject to change without notice.
6/2023 4157600g-en

Modification of this document is not permitted without written permission from Fluke Corporation.

